



GE Healthcare

Vivid S70N



M5Sc-D XDclear Active Matrix Single Crystal Phased Array Probe

Probe presets: Cardiac adult, pediatric, abdomen, fetal heart, transcranial, coronary, stress, LVO contrast

Biopsy guide: Multi-angle disposable with a reusable bracket



6S-D Phased Array Probe

Probe presets: Pediatric, cardiac, coronary, neonatal head, fetal heart, abdominal

Frequency: 2.4 - 8.0 MHz



12S-D Phased Array Probe

Probe presets: Pediatric, cardiac, coronary, neonatal head, rodent, vascular (incl. carotid, LEA, LEV, UE

Phased Array Probe

Frequency: 4.0 - 12.0 MHz



9L-D Linear Array Probe

Vascular (incl. carotid, LEA, LEV, UEA, UEV), musculoskeletal, nerves, small parts, thyroid
Biopsy guide: Multi-angle disposable with a reusable bracket

Linear Array Probe

Frequency: 2.4 - 10.0 MHz



C1-6-D XDclear-Single Crystal Curved Array Probe(Convex)

Probe presets: Abdomen, OB/GYN, urology, vascular, fetal heart

Biopsyguide: Multi-angle, disposable with a reusable bracket

Curved Array Probe (Convex)



P2D Pencil Probe

Probe presets: Cardiac

Frequency: 2.0 MHz



6Tc-RS TEE Probe

Probe presets: Cardiac, coronary, LVO contrast

Frequency Range: 3.0-8.0 MHz



6VT-D Active Matrix 4D Volume TEE Probe(4D and multi-plane only with 4D option)

Probe presets: Cardiac, LVO contrast, coronary

Frequency: 3.0 – 8.0 MHz



9T-RS TEE Probe

Probe presets: Pediatric

Frequency: 3.0 – 10.0 MHz